



Acoustic Rhinometer A1

Instructions for Use /
User Manual
Issue 4



USER MANUAL TABLE OF CONTENTS

A1 ACOUSTIC RHINOMETER SAFETY INFORMATION	PAGE 3
INTENDED USE	PAGE 5
DESCRIPTION	PAGE 5
STORAGE	PAGE 6
SPECIFICATIONS	PAGE 6
INSTALLATION	PAGE 11
INTRODUCTION	PAGE 13
HOW TO MAKE MEASUREMENTS	
ACOUSTIC RHINOMETER CALIBRATION	PAGE 14
PATIENT PREPARATION	PAGE 15
MAKING MEASUREMENTS	PAGE 17
FACTORS WHICH AFFECT ACCURACY	PAGE 17
ENVIRONMENTAL FACTORS	PAGE 18
A1 TROUBLESHOOTING	PAGE 19
INSTACAL	PAGE 21
A1TEST	PAGE 23
SPARE PARTS AND ACCESSORIES	PAGE 25
MAINTENANCE	PAGE 26
INSTRUMENT CLEANING AND DECONTAMINATION	PAGE 26
APPENDIX 1: GUIDANCE ON THE ELECTROMAGNETIC ENVIRONMENT IN WHICH TO OPERATE THE INSTRUMENT	PAGE 27
APPENDIX 2: TECHNICAL INFORMATION (RELATIVE to EN 60601-1:2006)	PAGE 31
Software Manual	PAGE 34

A1 ACOUSTIC RHINOMETER SAFETY INFORMATION



Read this Operating Manual before attempting to use the Instrument.

WARNINGS

This instrument is for indoor use only and it should be used as described in this manual.

The system must not be used in the presence of flammable gases or in an environment which is susceptible to explosions - beware of oxygen, dust, and anaesthetic gases.

To avoid risk of electric shock this equipment must be connected to a mains supply with a protective earth.

If your desktop PC or printer does not have a power supply approved for a patient environment, then an isolation transformer, which complies with BS EN 60601, should be used to power the PC, printer, and the A1 Acoustic Rhinometer. You must use the isolation transformer to ensure that the A1 Acoustic Rhinometer is in compliance with BS EN 60601.

As this is a mains-powered, PC-connected product it is advisable to not touch the patient while using the equipment.

The computer may not be connected to the internet.

The equipment should be positioned in such a way that it can be easily disconnected from the mains supply. The operation of the system can be safely terminated by switching off the device or removing the mains plug.

The applied part is the nosepiece, which is marked for single use. It should not be re-used, as no method of sterilisation can be guaranteed; disinfection materials may not be effective and may leave a residue, which could be harmful.

CAUTIONS

The nosepieces are made of a material, which may cause an irritation reaction in some patients. Use of the nosepiece should be discontinued in patients who exhibit such a reaction.

Nosepieces are a single-use item. Single-use items should not be re-used as they could transfer infections between subjects.

Federal (USA) law restricts this device to sale by or on the order of a physician.

The use of an A1 Acoustic Rhinometer near to sources of electromagnetic radiation, such as mobile phones, radio transmitters, x-ray equipment etc., may prevent it from functioning correctly.

The A1 Acoustic Rhinometer is a medical instrument, which has an electrical classification of Class I Type B and a Medical Device Directive classification of IIa. A Class I Type B device categorisation is used to describe an instrument which:

- a) Does not rely on basic insulation only to provide protection against electrical shock but is constructed in such a way that accessible metal parts cannot become live in the event of failure of the basic insulation, and
- b) Applied parts offer protection to the subject against electrical shock and in the event of a single fault condition arising, leakage current will be limited to less than 0.5 mA.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Contact GM Instruments at the address on page 10.

A1 should only be connected to other mains-powered devices such as computers and printers, which comply with EN 60950-1 and we also advise the use of an isolation transformer. Unless computers and printers built to EN 60950 are used, patient safety might be compromised.

Non-medical equipment such as computers and printers should be kept out of reach of subjects being tested as such equipment does not comply with medical safety standards. Refer to clause 16 of EN 60601-1:2006 to ensure compliance.

If the PC is allowed to go into sleep mode, the USB interface is powered down. When brought out of sleep mode, the PC does not re-initialise the USB interface, and therefore it is effectively not present. A solution is to remove the USB connection at the PC or instrument side then re-connect. That may be sufficient, but if not, save any results, then close down the instrument software and restart it.



Servicing can only be carried out by GMI approved and authorised personnel. No modifications are allowed.

Standardisation Document.

Customers are referred to the under noted International Standardisation Committee publication which contains recommendations on the use of Acoustic Rhinometers.

Acoustic Rhinometry: Recommendations for technical specifications and standard operating procedures.
Rhinology supplement number 10, December 2000 Ole Hilberg and Ole Pedersen

INTENDED USE

The A1 Acoustic Rhinometer is used to quantify nasal obstruction. The device is intended to be used by a Health Care Professional for non-invasive examinations of the nasal cavity using a sound pulse technology.

INDICATIONS FOR USE

The A1 is indicated for use where nasal obstructions are suspected, or when quantification of the nasal cavity before and after a procedure or challenge is required.

INTENDED USERS

The A1 Acoustic Rhinometer provides information to assist ENT Surgeons, Allergists, researchers into nasal function and those involved in sleep studies with information to guide their treatments and research. Typically, any medically trained technician, nurse, or doctor will be considered competent to use the device.

INTENDED POPULATION

The A1 Acoustic Rhinometer is intended to be used in paediatric and adult populations (including pregnant and breast-feeding women), regardless of gender.

BENEFITS

The A1 Acoustic Rhinometer has the following Benefits

1. Able to quantify the level of nasal obstruction
2. Provides a rapid non-invasive examination
3. Can be used for pre/post comparisons to assess outcome of intervention or therapy

DESCRIPTION of the A1 Acoustic Rhinometer

The A1 Acoustic Rhinometer allows a very rapid non invasive examination of the nasal cavity using a sound pulse technique.

The device comprises a hardware unit, with the sound tube attached by a cable, and software running on a computer. The sound tube contains a microphone, and is coupled to the test subject's nose with a nosepiece. The hardware unit is powered by a separate power supply, and connected to a computer via a USB cable.

Software running on the computer provides the user interface of the device, and sends and receives acoustic signals via the hardware unit.



Figure 1. The A1 Acoustic Rhinometer

STORAGE

The A1 Acoustic Rhinometer and its accessories should be stored within the following temperature and humidity ranges:

Temperature	-40°C to +60°C
Humidity	30% to 90% RH non-condensing

TECHNICAL SPECIFICATIONS

Only factory-trained personnel or engineers who are familiar with the standard EN 60601 can undertake servicing of the A1 Acoustic Rhinometer.





Circuit Diagrams will be made available to competent persons on request.

Performance (using standard nose)	Repeatability	Better than 2%
	Volume Accuracy (0cm to 5cm)	Better than 2%
	Minimum Area Accuracy (within 0cm to 5cm)	Better than 5%
	Area Range (max 20, default 10)	0.1 cm ² to 10 cm ²
	Distance Range (using standard sound tube)	0 cm to 12 cm
	Calibration	User Calibrated
Operation Environment	Temperature	+15 °C to +35 °C

	Relative Humidity	20% to 90% RH non-condensing
	Duty Cycle	Continuous
	Warm Up Time	5 minutes
	Supply	Universal voltage – external supply
	Power	10 Watts
Transportation and Storage	Temperature	-40 °C to +60 °C
	Relative Humidity	30% to 90% RH non-condensing
Mechanical Performance	Size	27 x 8 x 30 cm
	Weight	2 Kg

TABLE OF SYMBOLS USED










The following symbols appear on the A1 or the mains adaptor:





Symbol	Meaning	Socket Type	Location	Connected Part
	Refer to Instruction Manual ISO7010-M002	USB Connector Type B	Instrument Back Panel	Computer (Via USB port)
		5-Pin Din Socket	Instrument Back Panel	Mains AC/ DC Adaptor PCM50UT04*
	Type B Applied Parts IEC 60417- 5840	JIS C 5432-Compliant 5-Pin Socket	Instrument Back Panel	Sound Tube
	Protective Earth IEC 60417 - 5019	—	Instrument Back panel - Internal	—
	Direct Current IEC 60417 - 5931	—	Mains Adaptor AC/DC	—



For connected parts marked * only connect the accessories supplied with the instrument. These parts have been tested for use with the instrument for compliance to standards IEC 60601-1 and IEC 60601-1-2.

Symbols used on labelling and packaging:

SYMBOL	Indication
	Manufacturer
	EU Authorised Representative
	Date of Manufacture
	Catalogue Number
	Serial Number
	Medical Device
	Consult Instructions for Use
	CE Mark
	Do not re-use

	Non- sterile
	Temperature Limit
	Humidity Limitation
	Unique Device Identifier

ADDRESS AND CONTACT DETAILS:

MANUFACTURED BY:

GM Instruments Ltd, Block 1 Annickbank Innovation Campus, Annick Road, Irvine, KA11 4LF, UK.

TEL: +44 (0)1294 554664

EMAIL: sales@gm-instruments.com

WEBSITE: www.gm-instruments.com

AUTHORIZED REPRESENTATIVE:

Advena Ltd, Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013, Malta

INSTALLATION

INSTALLING THE A1 SOFTWARE

Installation of the A1 software is described fully in the accompanying manual titled:

A1 Acoustic Rhinometer Software Manual

INSTALLING THE A1 UNIT

The A1 Acoustic Rhinometer can operate on any PC which runs a Windows based operating system which has a free USB socket.

Connections between the A1 box, its sound tube, its power module and the PC are by means of polarised connectors which can only be inserted in the correct socket and with the correct orientation.

Connect the sound tube and the power module prior to powering up the A1, **but do not connect the USB cable until the A1 software and A/D board driver software are installed.**


N.B. Units which comply with EN60950 should be used.

HARDWARE IDENTIFICATION

Calibration plug----- supplied with the A1 Clinical and A1 Clinical/Research
Artificial nose and Straight Tube----- supplied with the A1 Clinical/Research



A1 rear view and connections



NOTE: LINK SOCKET ONLY USED WHEN CONNECTED WITH NR6 IN NARIS CONFIGURATION

INTRODUCTION

Hardware and software have been provided with the A1 Acoustic Rhinometer to allow for the determination of **Nasal Obstruction** by plotting Nasal Area against Distance using Acoustic Reflections.

The method is based on the following principle:

A sound pulse propagates in a tube and enters the nasal cavity through a nosepiece, where it is reflected by local changes in cross-sectional area. The incident and reflected signals are then measured by a microphone in the sound tube.

From these measurements, it is possible to calculate the cross-sectional area as a function of the distance into the cavity by use of algorithms developed by Ware and Aki (1969).

By integration of this curve the volume of the nasal cavity is calculated.

The method has been developed for use in the tracheo-bronchial system by Jackson et al. (1977), and for use in the nose by Hilberg et al. (1989).

For further information, the reader should consult these publications:

1. A.C. JACKSON, J.P. BUTLER, E.J. MILLET, F.G. HOPPIN Jr., and S.V. DAWSON.
Airway geometry by analysis of acoustic pulse response measurements.
J. Appl. Physiol. 43(3):523-536,1977.
2. J.A. WARE, AND K. AKI.
Continuous and discrete inverse scattering problems in a stratified elastic medium
I. Plane waves at normal incidence.
J. Acoust. Soc.Am.45:911- 921,1969.
3. O. HILBERG, A.C. JACKSON, D.L. SWIFT, and O.F. PEDERSEN.
Acoustic rhinometry, evaluation of nasal cavity geometry by acoustic reflection
J. Appl. Physiol. 66:295303,1989.

HOW TO MAKE MEASUREMENTS


The A1 can be used by any medically-trained technician, nurse or doctor, who has either read this manual or who has been trained by a competent authority in its use.

Having become familiar with the hardware and software you are now ready to record measurements.

CALIBRATION

Prior to use the A1 device **MUST** be calibrated for correct operation.

To perform a calibration, switch on the A1, put the calibration plug in the sound tube, click on the Acoustic icon on the main toolbar and then on the **Calibrate** button as shown below.



A1 Acoustic Rhinometer [USB-201] calibrated will show when completed.

A1 sound tube with calibration plug fitted, ready for calibration.

The system is now ready for use.



PATIENT PREPARATION

The patient should be prepared prior to testing, by having them wait in a quiet, stable environment for 15 to 20 minutes. If your interest is surgical and you want to look at structure, they should be decongested and asked to blow their nose prior to measurement.

If your interest is in monitoring response to a decongestant or to an allergen, again the patient should be asked to blow their nose prior to measurement.

The nosepieces supplied with the system include 15 conical and 5 pairs of each size of anatomical probes (medium & large). These are handed and colour-coded, i.e. one is for use on the right side (red) and one for the left side (blue).

N.B. The caution statement on page 3 at the front of this manual.



The choice of nosepiece depends on which will give the best seal with the least distortion.

The conical nosepieces are suitable for use with a nose which has a small round opening, while the anatomical nosepieces are designed to sit on the outside of a nose with a long narrow opening, and are shaped to make it easy, in most circumstances, to get a good seal. In addition, they have a flange which a soluble gel can be put onto to reduce any remaining gap.

MAKING MEASUREMENTS

- 1) Have the patient sit back in a chair (some workers advocate use of adjustable forehead and chin rests), select and fit a nosepiece, and apply the sound tube (perhaps suitably prepared with gel) to the side to be measured. The sound tube can be held by the patient, or by the clinician but we believe this is best left to the patient, with one hand holding the sound tube near the bottom and one near the top.
- 2) When suitably positioned, ask the patient to close their mouth and **gently** breathe in and out through their nose to check for leakage between the sound tube and nose.
- 3) If OK, ask them to open their mouth and return to mouth-breathing. When ready to make the measurement, ask the patient to take a breath in through their mouth and hold their breath. Click on the red **Acquire** button, or middle-click on the mouse with the cursor inside the capture window to sample the data. This process can be repeated until at least 3 consecutive results show little variation. If the patient can co-operate, the measurements may be made in rapid succession (with 600mS or so between samples) during a single breath-hold manoeuvre.
- 4) If the results are satisfactory, select the other side (L or R), and repeat the process, remembering to change the nosepiece, if anatomical.
- 5) The test results are automatically added to the patient's results file.
- 6) Click on **Print** to print out the currently highlighted curves, or alternatively click on **Report** to select the curves from the displayed list for printout. If you want a full page print out of a particular test, you can highlight it and click on the arrow on the right of the printer symbol and select **Print on full page**.
- 7) The results can be saved by clicking on **Save**, confirming the file name to save to, and then clicking on **OK**.

SWITCHING OFF

The software can be closed by clicking on **File** and then selecting **Exit** from the drop-down menu. The hardware can be powered down by unplugging, or by switching off the mains plug.

N.B. If you allow the PC to go to sleep or hibernate, the USB is disabled and will not re-initialise unless you either:

- a) Restart the PC or
- b) Close the A1 acquisition screen (if open), unplug the USB cable and then re-insert it.

We recommend that you disable the sleep/hibernate function, under Windows.

FACTORS WHICH AFFECT ACCURACY

Measurement accuracy and repeatability will depend on the following:

Calibration:	Perform a calibration at least twice a day
Temperature:	Use equipment in controlled conditions between +15°C & +35°C
External noise:	Use equipment in conditions where the background noise is below 65 dB
Angle of probe (relative to the head):	Aim for floor of nose
Rigidity of probe:	Doesn't move
Effect of probe on the nose:	Avoid distortion
Patient co-operation:	Hold breath during measurements
Nosepiece/nose seal:	Use gel and careful positioning to ensure a good seal.
Check for leaks/Distortion:	Repeat test at least twice if not 3 times and check for variation, i.e. apply/remove probe 2 or 3 times, testing each time and compare results.
Electromagnetic disturbances:	Exposure to excessive fields from nearby RF transmitters may the equipment not to function or measure correctly. If the measurement trace is not parallel to the horizontal axis in the first range of -4.5 to 0 centimetres, then the results should be discarded as invalid.

The Acoustic Standardisation Committee have now reported on the measurement and this is referred to on page 4 of this manual.

ENVIRONMENTAL FACTORS

Temperature Effects

An increase in the temperature of the gas through which the sound pulse travels lowers its density and therefore the rate at which it travels through the gas. The overall effect is in the range of 3% / 20°C. Distance is therefore slightly overestimated if there is an increase in temperature.

Altitude Effects

The effect of altitude is significant on the density of a gas and, therefore, on the speed of sound at an altitude of 1000m. As such distance will be overestimated by 7%.

Electromagnetic Disturbances

Exposure to excessive fields from nearby RF transmitters may the equipment not to function or measure correctly. If the measurement trace is not parallel to the horizontal axis in the first range of -4.5 to 0 centimetres, then the results should be discarded as invalid.

A1 TROUBLESHOOTING

If the PC is allowed to go into sleep mode the USB interface is powered down. When brought out of sleep mode the PC does not re-initialise the USB interface, which means that it is effectively not present. A solution for this is to unplug the USB connector at either the PC or instrument side and then re-connect. That may be sufficient, but if not, save any results, close down the instrument software, and restart it.

START UP

An instance of the application is already active

- The software is already open in another window

Error adding font GM.TTF or the file GM.TTF could not be loaded

- The font has not been added to your Windows font folder

License key is missing or invalid!

- Either no license key is present or it does not match the program

License key does not match product!

- The key file in the program folder does not match the program

FILE MANAGEMENT

The requested file cannot be loaded in this application

- The record file you tried to load is not matched to this program

File saving attempt failed!

- The record file has not been saved --- is the folder valid and do you hold the correct permissions to save to it?

Target file already exists!

- You are attempting to overwrite an existing file with this name

Directory "%s" does not exist

- The directory (folder) you are trying to save to or load from does not exist.

No filename given

- You have clicked on save without specifying a file name

Export to "%s" has failed!

- The data export you have set up has failed. Perhaps the folder name is wrong?

ACOUSTIC ACQUISITION

The hardware board with number %d is not recognized

- The USB A/D board you have specified is not installed (open InstaCal to register the board)

Error Code 10 - Invalid sampling rate (MCC UL Error 24 - Invalid sampling rate specified)

- The sampling rate specified in the program setup is too high for the USB or PC hardware to cope with - alter the basic sampling rate in the Devices tab (File/Settings/Devices)

Error Code 11 - The selected gain is not supported (MCC UL Error 30 - Invalid range specified)

- This applies to the setting up of A1Test.exe, where the gain parameter has been incorrectly entered

Error Code 12 - Sampling rate too high (MCC UL Error 29 - Overrun)

- As 10 above

Error Code 1001..1008 - Arithmetic error during computation

- May be a PC memory issue

Error Code 1009 - Signal below trigger level

- Either there is no click, it is weak, or the microphone or cable is faulty.
- (Use A1Test.exe to measure signal strength) If appropriate, VR2 may be increased if the signal is only a little below the trigger level or the trigger level setting can be reduced to below the measured signal level in the Acquisition (adv) tab. Details on how to use A1test.exe can be found below

Error Code 1050 - Error saving calibration data!

- The new calibration file generated has not been saved. Either the folder name specified during setup is not present, or you don't have permission to write to it

Error Code 2000 - Acquisition aborted

- Software/hardware issue has prevented the program from running

It is important that the software is loaded in the sequence suggested in the installation section of the **A1 Acoustic User Software Guide**.

If issues arise when trying to use the system, the following can be used to check that the system has been correctly set up, and to diagnose any problems:

InstaCal


Navigate to the MCC (Measurement Computing Company) folder using **My Computer** and double-click on the program **Inscal32.exe** - C:/Program Files(x86)/Measurement Computing/DAQ.

N.B. In a 32-bit system this will be C:/Program Files/Measurement Computing/DAQ as a 32-bit system won't have an x86 folder.

The opening dialog will show which boards are installed.

- The board you installed should be numbered 0 (or the corresponding number in the A1 software, if for some reason you have multiple devices installed)
- Select **Test** from the top bar and **Analog** from the drop-down menu, then **CH0**, signal source **External** and A/D range **±10V (for 201 USB interface)**

The line you see should lie on 0 Volts and should respond to talking into the open end of the sound tube – if this is ok then the microphone and connections between A1, the A/D card and the PC are ok.



EXPLANATION OF HOW InstaCal WORKS

Each time InstaCal is opened, it looks to see what Measurement Computing Company interfaces it can find. The first one it finds it allocates to **board 0** and the next to **board 1**, and so on. It records the board type and serial number and automatically edits a .ini file to include this information.

If you change one acoustic system for another, it will have the same board type, but it will have a different serial number, which means that the USB will not be recognised and the driver will not be loaded.

If you run InstaCal again, it will look for Measurement Computing Company interfaces. It will realise that the one already recorded is no longer present (by comparing serial numbers) and offer to remove it. When you answer yes, InstaCal will offer to record the new USB interface and when you click **OK** it will edit the .ini file to include the new board type and serial number. If board 0 is free, it will automatically allocate it to board 0, otherwise it will allocate it to the lowest available board number.

It is for these reasons that you can't only unplug one A1 instrument and replace it with another. You must run InstaCal to let it recognise that the first one is no longer present, delete it, and replace it with the new one.

If you have two USB devices connected, both can be registered by InstaCal one as board 0 and one as board 1. Provided you go into the software of one of them (**File -> Settings -> Device**) and change it to board 1 then you can easily switch between each system without having to run InstaCal.

A1TEST

In the A1 folder, you will find a program called **A1TEST.EXE** – this folder will be wherever you declared as the location during the installation process.


Double-clicking on the name will start the application.

For board type **USB-201** select Gain = ± 10 volts, Digital Port = AUX, Sample Rate = 100000.

Click on **Sample**. You should see the microphone response on the screen. If you point the mouse cursor under the trace, it will tell you the amplitude of the trace at each point.

The baseline should be between 1950 and 2050 units.

The value is only updated when you click on **Sample**. The picture below shows a typical display when the calibration plug is fitted to the sound tube and sample is selected.



The peak pulse size is also documented in the table on the left-hand side as **Maximum**.

It should peak at a level greater than that set within the A1 program for **Trigger**, and, using VR2, is normally set to **between 3200 and 3500**.

If the peak height is smaller than this the level can be increased by using the **gain control** resistor, VR2, on the A1 circuit board.


The trace above shows the microphone output plotted against sample numbers.

The first large peak shows the response when the click sound reaches the microphone on its way up the tube.

The second peak shows the reflected signal from the calibration plug returning back down the tube and reaching the microphone again.

Each pair of pulses that follow show the event repeating, as the sound travels up and down the closed tube.

If the calibration plug is removed, the trace will change as shown below and you will be able to hear the click when **Sample** is selected.



Selection of what constitutes incident and reflected waves is controlled by settings within the A1 program for:

incident wave window,
reflected wave start point and
reflected wavelength.

These are set for a particular length of sound tube and should not be altered.

SUPPLIED PARTS

INSTRUMENT AND ASSOCIATED PARTS	GMI CODE	QUANTITY
A1 Acoustic Rhinometer	A1-C	ONE
USB Cable	GM-USB	ONE
Sound Tube	A1-STM	ONE
Calibration Plug	A1-CP	ONE
External Power Supply	A1/PR	ONE
User and Software Manual	DOC-A1-UM-001	ONE
A1/MCC Software	A1-C-SW	ONE
*Large Anatomical Nosepieces	A1/NV1-PRL	5 Pairs
*Medium Anatomical Nosepieces	A1/NV1-PRM	5 pairs
*8 mm Conical Nosepieces	A1/NV1-PR8S	15 pcs
CLINICAL/RESEARCH VERSION		
Artificial Nose and Straight Tube	A1-AN.1	ONE

SPARE PARTS AND ACCESSORIES

NOSEPIECES:

*A1/NV1-PRM	Medium anatomically conformed nosepiece
*A1/NV1-PRL	Large anatomically conformed nosepiece
*A1/NV1-PR8S	Conical nosepiece

MISCELLANEOUS:

A1-CP	Replacement calibration plug
A1 Mic Sub	Replacement microphone and inner sound tube
GM-USB	USB Cable

If you have any questions about your Acoustic Rhinomanometer or require spare parts or consumables (part numbers listed above) then please contact your supplier or GM Instruments directly.

We will be able to advise you and give you help with any problem you may encounter.

Note: Parts listed above with (*) are supplied as non- sterile and due to their inertness have no expiry date.

MAINTENANCE MANUAL

OVERVIEW

The Acoustic Rhinometer should be calibrated each day using the built-in facility. In addition to this, the cables should be inspected regularly for signs of damage, in particular the region of the sound tube and its link to the PC. Do not use the equipment if damaged.

The nosepieces supplied are for single-use only. Care should be taken in use to prevent the passage of nasal secretions down the nosepiece and into the sound tube. Should this occur, the sound tube should be removed from the sound box, cleaned by one of the methods described below and then thoroughly dried prior to reconnecting it to the sound box.

N.B. Take care not to over tighten the fixing screws when re-fitting the tube as this could damage the tube.

INSTRUMENT CLEANING AND DECONTAMINATION

ENCLOSURE

Should the enclosure require cleaning for any reason, unplug it from the PC and wipe it with a damp cloth or cleaning wipes. Do not allow liquid to run into the enclosure.

SOUND TUBE

Under normal use, employing the disposable nosepieces supplied, no routine cleaning is believed necessary.

Should the sound tube enclosure require cleaning for any reason, unplug it from the A1 device and wipe it with a damp cloth or cleaning wipes. Do not allow liquid to run into the sound tube.

COPYRIGHT NOTICE

Copyright © G M Instruments Ltd, 2024. All rights reserved.
Copyright Protection.

Software and manuals are copyrighted, which means that your right to copy any part of them is limited by copyright law. Making copies or derivative works, except in certain limited instances and except for making copies for archival purposes or as an essential step in the utilisation of the computer program in conjunction with an instrument without prior authorisation is prohibited by copyright law and constitutes a punishable violation of such law.


Appendix 1: Guidance on the Electromagnetic environment in which to operate the instrument.


Guidance and manufacturer's declaration – electromagnetic emissions		
The A1 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1 should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Class B	The A1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The A1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC61000-3-3	Complies	

Appendix 1: Guidance on the Electromagnetic environment in which to operate the instrument (cont)

Guidance and manufacturer's declaration – electromagnetic immunity			
The A1 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -
Electrostatic discharge (ESD) IEC61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material,
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% U_T (>95 % dip in U_T) For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	<5% U_T (>95 % dip in U_T) For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the A1 requires continued operation during power mains interruptions, it is recommended that the A1 be powered from an uninterruptable power supply or a battery.
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Appendix 1: Guidance on the Electromagnetic environment in which to operate the instrument (cont)

Guidance and manufacturer's declaration – electromagnetic immunity			
The A1 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the A1, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance (<i>d</i>)</p> $d = 1.2\sqrt{P}$ <p><i>d</i> = 1.2√P 80 MHz to 800 MHz</p> $d = 2.3\sqrt{P}$ <p><i>d</i> = 2.3√P 800 MHz to 2.5 GHz</p> <p>Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in metres (m).</p> <p>Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
	6Vrms in ISM bands between 0.15-80MHz	6Vrms	
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

Proximity fields from RF wireless communications IEC61000-4-3	Frequency - Test Level 385MHz - 27 V/m 450MHz - 28 V/m 710MHz – 9V/m 745MHz – 9V/m 780MHz – 9V/m 810MHz - 28 V/m 870MHz - 28 V/m 930MHz - 28 V/m 1.72GHz - 28 V/m 1.845GHz - 28 V/m 1.97GHz - 28 V/m 2.45 GHz - 28 V/m 5.24GHz – 9V/m 5.50GHz – 9V/m 5.875GHz – 9V/m	Frequency - Test Level 385MHz - 27 V/m 450MHz - 28 V/m 710MHz – 9V/m 745MHz – 9V/m 780MHz – 9V/m 810MHz - 28 V/m 870MHz - 28 V/m 930MHz - 28 V/m 1.72GHz - 28 V/m 1.845GHz - 28 V/m 1.97GHz - 28 V/m 2.45 GHz - 28 V/m 5.24GHz – 9V/m 5.50GHz – 9V/m 5.875GHz – 9V/m	Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the NR6 is used exceeds the applicable RF compliance level above, the A1 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the A1</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

APPENDIX 2: TECHNICAL INFORMATION (relative to EN 60601-1:2006)

Warning: No modification to this equipment is allowed.

6.2 The A1 Acoustic Rhinometer is powered by an external universal mains power module. Power can be removed by unplugging the mains plug, or by switching power off at the mains socket.

It complies with the specification for Class 1 ME equipment, and is used along with a PC and printer, normally supplied by the customer.

N.B. The warning on page 3 about the PC and printer being compliant with EN 60950 and placing these items out of the reach of a subject

The applied parts are the nosepieces which are classified as a type B applied part.

6.3 The enclosure does not provide protection against the ingress of liquids and is therefore classified as IPX0.

6.4 There are no parts which are supplied sterile or which are required to be sterile.

6.5 The A1 Acoustic Rhinometer is **NOT** suitable for use in an oxygen rich environment.

6.6 The A1 Acoustic Rhinometer is rated for continuous use.

7.2.2 The A1 serial number can be found on the label on the back panel, and the software version number can be viewed on the start-up screen and by clicking on **Help** -> **About...** found on the top toolbar of the A1 software.

7.2.3 Information is given in the **Warnings and Cautions** section (page 3) of this manual, which the user must read prior to operating this instrument. Appropriate labels are fixed to the back panel.

7.2.4 Nosepiece packaging labels identify these as for single use only.

7.2.5 The A1 power supply module is marked PCM50UT04. No other external supply can be used, unless supplied by GM Instruments as a replacement.

7.2.6 The A1 power supply module provides DC voltages of ± 12 volts and + 5 volts from an input which can range from 100 Volts to 240 Volts AC at 50 Hz to 60 Hz. The power supply module is rated at 42W.

7.2.7 The A1 consumes 10W.

7.2.11 The A1 Acoustic Rhinometer is rated for continuous use.

7.2.17 Environmental conditions for transport and storage. Temperature: -40 °C to +60 °C
Humidity: 30 to 90% RH

7.9.1 The A1 can be used by any medically-trained technician, nurse or doctor, who has either read this manual or who has been trained by a competent authority in its use.

7.9.2.2 Warning: To avoid risk of electric shock, this equipment must be connected to a mains supply with protective earth.

7.9.2.3 The power module provided with the A1 (PCM50UT04) is considered to be part of the ME equipment.

7.9.2.5 The nose pieces are considered to be the applied part.

7.9.2.7 The equipment should be positioned so that it can be disconnected from the mains supply quickly and easily.

7.9.2.10 Error Messages --- see **A1 Troubleshooting** on page 19.

7.9.2.11 The A1 software can be closed by clicking on File -> Exit from the main toolbar. The A1 hardware can be switched off by removing the mains plug or by switching off the mains plug at the socket.

7.9.2.13 There are no parts which need routine maintenance or servicing other than general "housekeeping", such as keeping the equipment clean and cables free of twisting. Operational tasks such as calibration, and performance adjustments using the straight tube and artificial nose, are described elsewhere in this manual and in the sheets which accompany the tube and nose.

7.9.2.15 The applied parts (nosepieces) should be disposed of after use in line with your hospital or clinic's policy on disposal of potentially contaminated plastic parts. The A1 Acoustic Rhinometer can be returned to GM Instruments for dismantling and disposal, while your PC/Printer should be returned to your supplier for dismantling and disposal. These parts are subject to the WEEE directive and should not be disposed of in landfill.

7.9.3.2 The only parts which are interchangeable by service personnel are the power module and the mains cable. If these are changed, a PAT test to the current standard should be made prior to allowing release for use.

7.9.3.3 Circuit diagrams, component lists, and parts lists are available on request, along with email/telephone advice to service personnel trained on EN 60601-1:2006 and qualified to work on ME devices. Modification of the A1 Acoustic Rhinometer is not allowed.

7.9.3.4 If access to the A1 circuit board is required:

Disconnect the A1 from the mains supply by unplugging or switching off at socket

- Turn the unit over and remove the four feet on the base using a screwdriver

- Carefully turn the unit back over and remove the top panel
- The circuit board can now be accessed and, if voltage measurements are required, the mains power can be reconnected



Acoustic Rhinometer A1

User Software Manual and Installation Notes



USER SOFTWARE MANUAL TABLE OF CONTENTS

SOFTWARE INSTALLATION	PAGE 5
HOW TO CONFIGURE THE SOFTWARE	PAGE 39
STARTING A TEST	
TO PREPARE FOR MEASUREMENTS ON A NEW SUBJECT	PAGE 54
TO PREPARE FOR MEASUREMENTS ON A REPEAT SUBJECT	PAGE 54
TO MAKE A NEW ACOUSTIC RHINOMETRY RECORD	PAGE 55
TO COMPARE TEST RESULTS (L/R AND L/R-PRE AGAINST L/R-POST)	PAGE 57
TO DELETE A TEST RESULT	PAGE 60
BATCH TEST FACILITY (CLINICAL/RESEARCH VERSION ONLY)	PAGE 60
DATA EXPORT FACILITY (CLINICAL/RESEARCH VERSION ONLY)	PAGE 63

SOFTWARE INSTALLATION

COMPUTERS SUPPORTED

The A1 can be made to work on any PC which runs a Windows based operating system which has a free USB socket

N.B. Units which comply with EN 60950 should be used.

PRINTERS SUPPORTED

As these instruments operate in a Windows environment, print capability depends on you having installed a printer under Windows. Virtually any printer, which works under the version of Windows you have, will be suitable.

INSTALLATION SEQUENCE

The software comes on 2 CDs or 1 Thumb Drive (containing 2 folders). It is advisable to close any other programs you may be running while performing the installation as you may need to restart the PC during this operation.

INTERFACE BOARD SOFTWARE (Measurement Computing) Ensure the USB cable is NOT connected to the A1 unit.


If you have an MCC DAQ disk, put it in your drive and allow it to auto run, or open it and click on Install. You only need to install InstaCal/Universal Library software from this disk. Other options such as Tracer and manuals are not required.

If you have a Thumb Drive put it in a free USB socket and navigate to it. Open the MC folder, and click on icalsetup.exe. Accept the defaults offered.

The program InstaCal and other driver software will be loaded to a folder called Measurement Computing.

INTERFACE BOARD HARDWARE


The next installation operation is to have your PC identify the USB A/D converter, which is built into the A1 unit. To do this, connect the A1 to your PC using the cable supplied. (It is not necessary to power up the A1 or at this stage or load its software) Click on the **Start button**, scroll down to **Measurement Computing**, select it and scroll down to **InstaCal** and click on it. **InstaCal** will show no boards installed as shown below.




Wait for a few seconds as the USB interface is recognised and the drivers are loaded.

If this doesn't happen, click on **Install**, then **Refresh**. Once completed, the InstaCal screen will offer the opportunity to accept your USB board.

N.B. the serial number for the USB in each A1 is unique. The most recent units are for type 201, while older A1 units will have a 1208 USB module installed.



Click on **OK**



When InstaCal shows the screen above, (shown for a 201 module) USB installation is complete and InstaCal can be closed.


GMI SOFTWARE

Put the GMI CD in your drive or open the folder A1 on your Thumb Drive. Navigate to the file Setup-Acoustic_V3_X_X_XXXX.exe and click on it. Accept the defaults offered at each stage, until you are asked for the desired location for the .ini file – you will be offered the choice of the shared folder or the application folder. Unless you plan on using the system on a network, the application folder should be selected.

The A1 Acoustic Rhinometer program can now be started by clicking on:

The **Start** button and from the program list select Acoustic, then Acoustic.

The program opens a blank patient file as shown:



If you wish to change the fields listed in the Patient Record screen, this should be done at the beginning, before starting to save records. To do this, click on File, click on Close, click on File, click on Settings, select the Patient Record Tab and make the alterations you want, as described in the following pages.

If you wish to make measurements on a New subject, click on File, click on Close, then click on New.

If you wish to bring back a record of a Patient previously tested, click on Open.

HOW TO CONFIGURE THE SOFTWARE

Several elements within the program can be configured to suit your application.


N.B. Although you may wish to reconfigure these elements again later it should be noted that fields within the Patient record screen cannot be changed after storing results for a particular subject as the stored information may become inaccessible – see page 10 for more details.

All changes are made by clicking on **File -> Settings**.

N.B. The values shown on the following screens (and particularly in the acquisition screens) are examples only, as every system is individually configured.

A new window named **Settings** is then displayed which shows two super-tabs, each of which contains subtabs with user-configurable items. The two tabs are labelled **General** and **Acoustic**.

GENERAL > USER INTERFACE TAB




The **General > User Interface** tab allows you to:

- Show captions on the toolbar buttons
- Show the PC system time in the status bar
- Enable the middle mouse button to open the acquisition window

- Make the system adjust the graphics display following a test to allow space for the table of results appropriate to that test
- Make the system match the language set on the PC to any of the options supplied with the A1 software, when possible

GENERAL > PATIENT RECORD TAB



The **General > Patient record** information window can contain many fields. You can include as many or as few as you want, put them in any order, define what kind of data is acceptable for entry in any field and specify the size of each field.

The Patient Record tab contains the following elements:

List of fields is a window which shows the current labels already selected for inclusion. If you click on one of these it can be edited or deleted or its position moved within the group. When highlighted, the information relating to its structure is shown on the right-hand side of the window under the following headings:


Caption - the current label selected for editing

Maximum Length - the number of character spaces allocated in the **Patient Information** window against that name

Input mask (format) - the conditions applied to character entry in the associated **Caption** field

In the above example, **ID no.** is highlighted in the **Lists of fields** window. The label **ID no.** is shown in the **Caption** field and it can be seen that currently space is allowed for 15 characters. The applied mask was one which only allowed valid filename

characters, and no spaces or other characters. Masks can be selected or changed by clicking on the button on the extreme right of the **Input mask (format)** field.



The **Input masks** window allows selection of the conditions you want to apply and has space below for you to test the restrictions. If you wish to keep the currently selected field, click on **OK**, and then on **Apply**.

If you wish to discard any field, highlight that field and click on **Delete**. If you want to add a new one, click on **New**.

If you want to change the order in which they appear, highlight the one you want to move and then use the **Up** and **Down** arrow buttons to the right of the **Delete** button.

You can place more than one item on a line by placing it below the item you want it lined up with and clicking on **Position on same line as previous field**.


You can also hide a field from a printout by highlighting it and clicking on **Hide on printouts**.

A field can be used for the test file name by highlighting it and then clicking on **Use for suggestion of filename**.

N.B. The software will open in a blank patient record. If you wish to add or remove any fields from the Patient record screen or make changes to them you must first close the current record by selecting **File -> Close** from the main toolbar.

(This software enhancement has been made to stop patient files being accidentally corrupted – a side effect of a patient record being changed while open within the software).

GENERAL > STYLES TAB



The **General > Styles** tab allows you to set up colours and line styles for either the screen or the printout.

Target - here you can select **Application user interface (screen)** or **Report (printer)** for alteration

Plot Item - scrolling window which contains a list of the screen or printout elements whose colour, line style, or font can be changed

Color - colour palette, from which you choose the colour you want applied to the element selected in the item window

Style - displays the colour, line style or font currently set for the element shown in the item window

Pen Style - allows you to set different types of line for the selected element

Thickness - allows you to apply different line thickness to any line drawing element selected from the item window

Preview - window which lets you see the effect of a change


To make a change to the default settings select **Application user interface (screen)** or **Report (printer)** from the top dropdown menu, select an **Item** from the list, select your desired **Pen** colour and, if appropriate, **Pen Style** or **Thickness**. You can immediately see the effect of that change in the **Preview** window.

If you click on **Apply**, the change will be recorded and you can move to another element. If you click on **OK** instead of **Apply**, the change will be made, and recorded, but you will close the **Setup** window and return to the main screen.

Application User Interface (screen) & Report (printer)

Item Name	Item Description
Acoustic Plot Background (Testing screen only)	The background colour of the graphs viewed on the screen during testing
Axis Labels	Numbers marking the axes on the A1 graphs
Matrix	The non-zero grid lines for each graph
Curve 1 – Curve 10	Each of the captured curves shown on a graph at any given time
Acoustic MeanBase (L+R)	Colour of the left and right baseline curves when performing a comparison
Acoustic MeanComp (L+R)	Colour of the left and right comparison curves
Acoustic Goldnose	Curve relating to the Goldnose Acoustic Rhinometer calibration device
Graph Border / Graph Area	Colour of graph borders and background areas
Axes	Zero lines on the graphs
Results Values Text	All text in the tables
Subgraph Caption	Text in the “Remarks” box in the bottom right corner of the screen
Legend: Item Text/Box Frame/ Box Background	Colour of the frame, background within the frame and the text within the frame of a legend describing the information contained within a graph
Patient Record Background	Background of the main screen where patient details are entered


GENERAL > REPORT TAB




The **General > Report Tab** allows for the setting of a title which will appear on printouts and the alteration of printer margins. More detailed control of printer options is available in the normal Windows print driver dialogue box.

ACOUSTIC > DEVICES TAB

The **Acoustic > Devices** screen shown below displays the USB device number assigned to the instrument [0]. If you then click on the device (**A1 Acoustic Rhinometer**, in the below example) the buttons at the bottom of the window become active.




If you click on **Configure**, you will see:



When the USB interface is installed (**USB-201**), the window will show: the board number [0], the channel used for the microphone (which, by default, is channel 0), the selected gain setting (**±10 Volts**), the basic sampling frequency (100kHz), the upsampling option will be disabled, and the pause between each measurement (**700 mS**). These are the default settings for a **USB-201** interface and should not normally be changed.


N.B. Earlier versions of the A1 may have also used a USB-1208FS interface, in this instance the window will show: the board number [0] and the channel used for the microphone (0), the selected gain setting (**± 1.25 Volts or ± 10 Volts**), also displayed are the basic sampling frequency (**46kHz**), the data processing upsampling equivalent (**100kHz**) and the pause between each measurement (**700 ms**).



When you click on **Calibrate** in the test acquisition window, you should have the A1 powered up, with the calibration plug fully inserted into the A1 sound tube. An updated calibration file will be generated and saved to the default folder – the same directory to which the program was installed.

The **Add** and **Remove** buttons allow the interface to be changed, if, for instance, you have an NV1 Rhinospirometer, which you want to install on the same PC. Each interface will be assigned a different number by InstaCal, and, therefore, as the default setting for both instruments is board 0, one or other must be changed to allow both to work. Normally there is no need to alter these settings.


ACOUSTIC > ACQUISITION TAB



Side - sets the side of the nose which is, by default, to be tested first.

When acquiring acoustic data, the information can be obtained from a single measurement or from the average of a series of curves. The item labelled **Number of traces** specifies what is to constitute a single test.

ACOUSTIC > ACQUISITION (ADV) TAB



All the parameters set up within this dialog box are pre-set by GM Instruments, to give best possible performance on your system, and are stored on your program CD.

Every sound tube is different and, therefore, may have different setup parameters associated with it.

The **Distance shift** factor allows a small lateral movement of the measured data curve to give the best fit when comparing a measurement made with an artificial nose against the theoretical line produced by the profile of an artificial nose. It effectively compensates for variations in the speed of sound under different conditions.

Sound tube area, **Zero adjust**, **Stability**, **Starting point shift** and **Attenuation** factors are adjusted to provide optimum results when setting up the system using with the artificial nose and straight tube, as recommended by the standardisation committee.

The **Max distance** factor is normally set to **25 mS** to limit the size of the files which are created and saved. Normally, scaling of the horizontal axis is set to run from –5 cm to +15 cm, so the minimum requirement is a maximum distance of 20. If you want to measure further back, say to 25 cm, then the horizontal scale would need to be set to run from –5 to +25 and the maximum distance set to at least 30.

If an improper data message is produced during acquisition the **Trigger level** factor can be reduced to allow measurements to be made.

Data windows parameters are based on the physical dimensions of every sound tube and are, therefore, constant for every sound tube.

Nosepiece Length allows you to tell the software what nosepiece you are using to let it adjust the graphical display accordingly – this will normally be 4.5 cm.

ACOUSTIC > RESULTS TAB

Settings

- General
 - User interface
 - Patient record
 - Styles
 - Report
- Acoustic
 - Devices
 - Acquisition
 - Acquisition (adv)
 - Results
 - Plot
 - Reference data

Acoustic > Results

Volume calculation

Start distance: 0 cm

Stop distance: 5 cm

Comparison ranges: 0-5, 2-5 cm
Maximal 5 ranges: start-end. Separate ranges with comma (e.g. 0-5, 2-5)

Cross-sectional area minima

Start distance: 1 cm

Stop distance: 5 cm

Hysteresis: 0.015 cm²

Comparison distances: 2,4 cm
Maximal 5 values. Separate values with comma (e.g. 3, 5, 7)

Miscellaneous

☒ Include totals

☒ Include Minimal Cross-sectional Area (MCA)

OK Cancel Apply

Volume calculation start distance and **stop distance** specify the section of the nose over which you want to calculate a volume.

Cross-sectional area minima specifies the section of the nose over which you want to look for turning points, and **Hysteresis** specifies what constitutes a turning point.

Comparison ranges allows you to specify the volumes reported between up to 5 pairs of distances, when a comparison display is on screen.


Comparison distances allows you to specify up to 5 distances at which cross-sectional areas will be reported when a comparison display is on screen.

Miscellaneous

The **Include totals** tick box is offered here to allow the standard result screen, when showing right and left side volumes etc., to also give a total nose figure. This summation may not be appropriate if you want to compare two right side or two left side measurements, and can therefore be switched off here.

Although the system will automatically report MCA1 and MCA2 (the first and second Minimum Cross-Sectional Areas found, starting from the front of the nose), it can also automatically determine which is the smallest and report it as MCA, if the box is ticked.

ACOUSTIC > PLOT TAB




The scaling for the distance and cross-sectional area axes can be set up here, for the **Linear plot** and the **Logarithmic plot**, comparison and analysis view screens.

N.B. A logarithmic scale should not be assigned a minimum value of zero.

Miscellaneous

A tick box is offered to allow the side used to display the left/right results to be switched for comparison tests.

ACOUSTIC > REFERENCE DATA TAB



Validation of system performance can be made by comparing a theoretical **Standard nose** curve, with that obtained when making a measurement on a **Standard nose** model.

The A1 GMNose allows such a system check to be made. The **Acoustic > Reference data** tab shown above, allows the appropriate theoretical curve to be drawn on the screen in a colour and line style specified in the **General > Styles** tab, under control of check boxes on this screen.

The **Datafile** used under normal circumstances is called GMNOSE.DAT and can be selected from file by clicking on the button to the right of the **Datafile** window. The selected data file can be shown on screen if **Visible on screen** is selected and can be seen on printouts if **Visible on printouts** is selected.

The **Datafile** is drawn behind measured curves provided the 3rd tick box is selected.


STARTING A TEST

TO PREPARE FOR MEASUREMENT ON A NEW SUBJECT

1. All current patient information and test results are cleared when the program is started, or if, having performed other measurements, you click on the **New** button.
2. In either event the cursor is set to the first field in the patient record window.
3. You can then enter the required patient details. Pressing the tab key after entering information in each field will take the cursor to the next one.
4. To make measurements (acquire), click on the Acoustic icon on the top toolbar.

TO PREPARE FOR MEASUREMENT ON A REPEAT SUBJECT

1. Click **Open** and then, using the dialog box offered, go to the folder you store your results in and click on the appropriate subject's file name.
2. You will then be able to examine the patient details and check that you have the correct subject by clicking on the **Patient record** tab. Clicking on the **Tests** tab will show you all previous test dates and results for that patient.




TO MAKE A NEW ACOUSTIC RHINOMETRY RECORD

Click on the Acoustic icon on the top bar and one of the following two windows will appear:


- 1) If no calibration file is found you will be prompted to **Calibrate**. With the A1 switched on and the calibration plug securely fitted to the end of the sound tube, click on the button marked **Calibrate**.



- 2) If a calibration file is found, you will see the screen below and can proceed to testing, or alternatively you can refresh the calibration by ensuring the A1 is switched on and that the calibration plug is correctly fitted, and clicking on the button marked **Calibrate**.



- 3) On the top bar, the number of **traces** representing a test can be changed from the default number (4). If you have a clinical/research version you can also switch on **Batch mode** to allow result validation (see page 28 for more details). The facility **Trigger on zero flow** is a specialised feature used in conjunction with a research-based combined acoustic and rhinomanometer sound tube.
- 4) Confirm that the side you want to test is selected (**left or right**), and once the patient is connected, click on the button with the red spot on it or, provided the mouse cursor is within the acquisition screen and you have this setting selected, your mouse middle button to begin the acquisition.
- 5) Once the pre-set number of measurements has been taken the system will stop taking in data, will show you the resultant curves, each of which is drawn in a different colour, and will add the test to the test list for that subject.




- 6) If any curve is obviously bad it can be de-selected by clicking in the appropriate check box in the **Traces** section on the right of the screen.
- 7) Comments can be added in the box found on the bottom right hand side of the screen by clicking **Edit...**, or removed by clicking on **Clear**.

- 8) If required, additional testing can be done, or the current test can be printed by clicking on **Print**, or it can be saved by clicking on **Save** – if a full screen printout is required, the desired test should be highlighted and **Print on full page** should be selected by clicking on the arrow next to the **Print** icon.
- 9) The **Output** tab can be used to create a printout consisting of a number of curves selected by clicking in the check box on the left of the test list.
- 10) Additional facilities available from the screen includes:
 - Import** --- This allows you to merge data from two saved files
 - Export** --- This allows the raw data to be sent to a text file
 - Copy** --- This allows the calculated values to be put on the clipboard, ready to be pasted into a spreadsheet

TO COMPARE TEST RESULTS (L/R AND L/R-PRE AGAINST L/R-POST)

If you want to create a complete nose record from Left-side and Right-side results:

- 1) Select the **Tests** tab
- 2) Click on a Left test using the left mouse button
- 3) Hold down the ctrl key and click on the desired Right test using the left mouse button - both tests will now be highlighted. Release the ctrl button, right-click on either test result and click on **Create comparison test**
- 4) A new test entry is created which shows the Left curve and the Right curve, which can be displayed in a number of different ways by selecting different toggle buttons, which are found below the graph



The above illustration shows a L/R comparison test, displayed in **Report** style.


If you want to create a test which compares two Left side or Right-side results:

- 1) Select the **Tests** tab
- 2) Click on the test you want to use as the baseline test using the left mouse button
- 3) Hold down the ctrl key and click on the test you want to be the comparison test using the left mouse button. Both tests will now be highlighted. Click on one of the tests using the right mouse button and from the drop-down menu, click on **Create comparison test**
- 4) A new test entry is created called **COMP**, which shows both curves. This can be displayed in a number of ways, which can be toggled via buttons beneath the central window

It is also possible to compare a pair of total nose records. In this case one pair (L+R) record is designated as the baseline result and a second pair (L+R) is designated as a comparison result (colours and line styles for these are specified in the **Style** tab).

Left-click on one L+R curve to designate it as the baseline, then hold down the ctrl key and click on the L+R curve you wish to compare against. Right-click on either test and select **Compare** to produce a new record entry.

This new record is called an **ANALYSIS** record. When selected, a detailed examination of the two sets of results is possible, and the % change at different points is calculated.



TO DELETE A TEST RESULT


Select the test result to be deleted from the test list using the left mouse button to view it on screen, then right-click on the test and, from the menu which appears, select delete.

BATCH TEST FACILITY (CLINICAL/RESEARCH VERSION ONLY)

Validation of test results is a facility which has proved to be invaluable in many situations. In acoustic rhinometry result variation is normally caused by distortion of the nose or leakage in the connections to it.

The batch test facility has been incorporated into the Clinical/Research version of the program in order to allow users to check for these sources of result variation.

The program offers the possibility of switching on the **Batch** process when you go into the acquire screen, by clicking on **Batch mode** in the top right corner of the window. An additional section of window is then added to the acquire screen, where a number of pieces of information will be added as test results are taken.



At the end of the first test **Start batch** will be shown in the progress dialog box. At this point the subject should be disconnected from the equipment, and then reconnected.

The red button in the bottom left corner should be clicked on again to start part two. Once completed, a coefficient of variance figure (**CV**) will be shown in the progress box, while additional information such as the mean resistance, standard deviation between values, and the **CV** figure are shown in the results section.

If the **CV** figure is low then there is a strong likelihood of the values being accurate, as the chance of creating the same distortion or leakage twice in succession is remote.

However, additional tests can be added to the batch (with the subject disconnected between each) if you want further re-assurance.

The batch values can be printed directly without any possibility of modifying the data by clicking on the batch print button.

When you close the batch test window, the tests are transferred to the normal test list for further examination and permanent storage.

DATA EXPORT FACILITY (CLINICAL/RESEARCH VERSION ONLY)

STANDARD PRINTOUTS

Graphical and numerical information can be outputted in the form of printed records either directly from the **Tests** screen by selecting the records required and clicking on the **Print** button, or as a composed report by clicking on the **Output** tab and selecting the records you want on your report. These can be viewed on the right of the screen prior to clicking **Print**.

SCREEN TRANSFERS

A single tap on your keyboard's **Print Scrn** key will transfer the contents of the screen to the clipboard, which can then be pasted into other record-keeping systems.


DATA OUTPUT

Numerical information can also be outputted from Clinical/Research versions of A1 in two ways.

- 1) Cross-sectional area values at pre-set time intervals:


This is done by saving a record or group of records to a file of type .DAT.

The .DAT format saves cross-sectional area values sequentially at a fixed distance interval, which depends on the acquisition rate. At the standard sampling rate of 100kHz, area values are presented every 1.73 mm along the distance axis.



Click on **Save**, select files of type .DAT, enter the file name, if not entered automatically, and then click on **Save**. The saved file can then be viewed in Notepad or put into other software for processing.

- 2) The Clinical/Research version of the software includes an export facility, which can be found under the **File** drop-down menu, which, when clicked, presents you with a setup dialog box:



This allows you to select what patient information you want to export to a file of type .EXP, and whether the information should be calculated values, raw data, or a mixture of both.

As a companion to the .EXP file, a .DEF file can be created at the same time which names each item of data in the .EXP file.

Both types of file can be opened within commercial packages such as Microsoft Excel and subsequently processed in any way you wish.

Approved version: 4

Title: Instructions for Use - A1 Acoustic Rhinometer

Approved By:

General Manager - Joe O'Neill

Thu Aug 7 17:33:36 GMT 2025

Approved and Signed

m48ZDzw6bf4ILRtim06pUvhnnek

Quality Engineer - Rachel Fraser

Fri Aug 8 10:36:07 GMT 2025

Approved and Signed

NNI19XLJej8GSA1RNSdXyKEU1mY

Approval History:

Version	Approved on	Status	Issued by
4	Fri Aug 8 2025	Approved	Rachel Fraser
3	Tue Jul 15 2025	Superseded	Rachel Fraser
Version	Approved on	Status	Issued by