



PEAK NASAL INSPIRATORY FLOW (PNIF) METER



Product Information

The GM Instruments PNIF is a simple, portable inspiratory flow meter, which measures PNIF (peak nasal inspiratory flow) through the nasal passage; providing information of nasal obstruction. Objective measurement of nasal obstruction is useful in the diagnosis of a number of conditions and establishing efficacy of nasal decongestants. Considered as an inexpensive alternative

- Measures inspiratory flow through the nose
- Accurate assessment of nasal congestion
- Individually calibrated
- Can be used during assessment of allergic rhinitis

Technical Specification:

Dimensions	158 x 44mm - w/o mouthpiece 299 x 138 x 106mm - with mask
Weight	131g (with mask)
Measurement Range	30 - 370 L/min
Accuracy	10% or ± 10 L/min
Repeatability	5 L/min

Product Codes:

PNIF with Medium mask	PNIF
Single use Child Mask for PNIF Meter	PNIF-SC
Single use Adult Mask for PNIF Meter	PNIF-SA
Reusable Child Mask for PNIF Meter	PNIF-RC
Reusable Adult Mask for PNIF Meter	PNIF-RA